Technical Data Sheet

InVivoMAb anti-mouse MHC Class II (I-A/I-E)



Lot Specific Information

 Lot Number:
 Lot Specific*

 Volume:
 Lot Specific*

Concentration: Lot Specific* (generally 4 to 11 mg/ml) *

Total Protein: Lot Specific*

*This information will be noted on the certificate of analysis that ships with this product.

Product Information

Reported Applications:

Catalog Number:BE0108Clone:M5/114Isotype:Rat IgG2b

Recommended Isotype Control(s): InVivoMAb rat IgG2b isotype control, anti-keyhole limpet hemocyanin

Recommended Dilution Buffer:InVivoPure pH 7.0 Dilution Buffer
Immunogen:

Activated C57BL/6 mouse spleen cells

in vivo MHC II blockade

Functional assays Immunofluorescence Western blot Immunoprecipitation Flow cytometry

PBS, pH 7.0

Formulation: Contains no stabilizers or preservatives

<2EU/mg (<0.002EU/µg)

Determined by LAL gel clotting assay

Purity: >95%
Determined by SDS-PAGE

Sterility: 0.2 µM filtered

Production: Purified from tissue culture supernatant in an animal free facility

Purification:Protein GRRID:AB_10949298Molecular Weight:150 kDa

Description

Endotoxin:

The M5/114 monoclonal antibody reacts with mouse MHC Class II haplotypes I-A^b, I-A^d, I-A^d, I-E^d, and I-E^k. The antibody does not react with I-A^f, I-A^k, or I-A^s haplotypes. The M5/114 antibody is reported to inhibit I-A-restricted T cell responses.

Shelf-life and Storage

Store at the stock concentration at 4°C . Do not freeze.

All Bio X Cell antibodies have a guaranteed shelf-life of one year from the date of customer receipt when stored as recommended. It is not uncommon for a floccule or precipitate to appear during storage. The floccule is typically buffer salts precipitating out of solution or a small bit of protein aggregation. For information on how to remove floccules or precipitates see our FAQ's at bxcell.com/faqs.

Protocol Information

Since applications vary, each investigator should use the application references as a guide to help estimate the appropriate dose or concentration. The dose or concentration can be further optimized experimentally in a dose response or titration experiment.

Application References

For a complete list of references, visit https://bxcell.com/product/mhc-class-ii/#references or scan the QR code below.

Bio X Cell, Inc.

bxcell.com

1.866.787.3444

customerservice@bxcell.com

Conditions: For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

Bio X Cell, Bio X Cell Logo and all other trademarks are the property of Bio X Cell, Inc. © 2020 Bio X Cell



Bio X Cell, Inc.

bxcell.com 1.866.787.3444

customerservice@bxcell.com

Conditions: For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.