## **Technical Data Sheet**

# InVivoPlus anti-mouse IFNy



Attention: Use of this product constitutes an agreement to Bio X Cell's Terms and Conditions which are included with this product in print and can also be found at <a href="https://bxcell.com/terms-and-conditions/">https://bxcell.com/terms-and-conditions/</a>.

## **Lot Specific Information**

Lot Number:Lot Specific\*Volume:Lot Specific\*

Concentration: Lot Specific\* (generally 4 to 11 mg/ml) \*

Total Protein: Lot Specific\*

\*This information will be noted on the certificate of analysis that ships with this product.

## **Product Information**

Catalog Number: BP0055 Clone: XMG1.2 Isotype: Rat  $\lg G1$ ,  $\kappa$ 

**Recommended Isotype Control(s):** InVivoPlus rat IgG1 isotype control, anti-horseradish peroxidase

Recommended Dilution Buffer: InVivoPure pH 8.0T Dilution Buffer

Immunogen: Recombinant mouse IFNy

*in vivo* **IFNγ neutralization** *in vitro* **IFNγ neutralization** 

Reported Applications: ELISPOT

Flow cytometry Western blot

Formulation: PBS + 0.01% Tween, pH 8.0

Contains no stabilizers or preservatives

Endotoxin: <1EU/mg (<0.001EU/µg)

Determined by LAL gel clotting assay

Purity: >95%

Determined by SDS-PAGE

Sterility: 0.2 µM filtered

**Production:** Purified from tissue culture supernatant in an animal free facility

Purification: Protein G

<5%

Aggregation: Control of the part of the pa

# **Murine Pathogen Test Results**

Mouse Norovirus: Negative, Mouse Parvovirus: Negative, Mouse Minute Virus: Negative, Mouse Hepatitis Virus: Negative, Reovirus Screen: Negative, Lymphocytic Choriomeningitis virus: Negative, Lactate Dehydrogenase-Elevating Virus: Negative, Mouse Rotavirus: Negative, Theiler's Murine Encephalomyelitis: Negative, Ectromelia/Mousepox Virus: Negative, Hantavirus: Negative, Polyoma Virus: Negative, Mouse Adenovirus: Negative, Sendai Virus: Negative, Mycoplasma Pulmonis: Negative, Pneumonia Virus of Mice: Negative, Mouse Cytomegalovirus: Negative, K Virus: Negative,

#### **Description**

The XMG1.2 monoclonal antibody reacts with mouse IFNγ (interferon gamma) a 20 kDa soluble pleiotropic cytokine and the sole member of the type II class of interferons. IFNγ is primarily produced by activated lymphocytes including T, B, NK cells, and ILCs. IFNγ exerts immunoregulatory, anti-proliferative, anti-viral, and proinflammatory activities and plays an important role in activation, growth, and differentiation of T and B lymphocytes, macrophages, NK cells and other non-hematopoietic cell types. Additionally, IFNγ induces the production of cytokines, Fc receptor, and adhesion molecules and up-regulates MHC class I and II antigen expression by antigen presenting cells during an immune response. IFNγ has also been shown to modulate macrophage effector functions, influence isotype switching and induce the secretion of immunoglobulins by B cells. IFNγ signals through the IFN gamma receptor which exists as a heterodimer composed of CD119 (IFN gamma receptor 1) and AF-1 (IFN gamma receptor 2). The IFNγ receptor is expressed ubiquitously on almost all cell types with the exception of mature erythrocytes. The XMG1.2 antibody is a neutralizing antibody.

## **Shelf-life and Storage**

Bio X Cell, Inc.

bxcell.com 1.866.787.3444

customerservice@bxcell.com

customerservice@bxcell.com

Conditions: For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

Bio X Cell, Bio X Cell Logo and all other trademarks are the property of Bio X Cell, Inc. © 2020 Bio X Cell

Store at the stock concentration at 4°C. Do not freeze.

All Bio X Cell antibodies have a guaranteed shelf-life of one year from the date of customer receipt when stored as recommended. It is not uncommon for a floccule or precipitate to appear during storage. The floccule is typically buffer salts precipitating out of solution or a small bit of protein aggregation. For information on how to remove floccules or precipitates see our FAQ's at <a href="mailto:bxcell.com/faqs">bxcell.com/faqs</a>.

## **Protocol Information**

Since applications vary, each investigator should use the application references as a guide to help estimate the appropriate dose or concentration. The dose or concentration can be further optimized experimentally in a dose response or titration experiment.

## **Application References**

For a complete list of references, visit <a href="https://bxcell.com/product/invivoplus-anti-m-ifngamma/#references">https://bxcell.com/product/invivoplus-anti-m-ifngamma/#references</a> or scan the QR code below.



Bio X Cell, Inc.

bxcell.com

1.866.787.3444

## customerservice@bxcell.com