




Bio X Cell

Your Trusted Supplier of *in vivo* MAbs

bioxcell.com



Your Trusted Supplier of *in vivo* MAbs

Accelerate your research from pilot to pre-clinical phases

Rely on us to deliver the largest catalog of monoclonal antibodies formulated for *in vivo* use and provide customizable antibody production services - backed by over 25 years of experience.

"When anyone asks me for CTLA-4 antibodies, I always refer them to you. I look forward to another 20 years."

- James P. Allison -
Corecipient of the 2018 Nobel Prize in Physiology or Medicine

Discover the Bio X Cell Advantages

Extensive Portfolio: Over 500 unique antibodies that are ultra-pure, low in endotoxin, and preservative-free.

Extensive Reference Library: Leverage over 20,000 product citations to craft your experiments more effectively.

Exceptional Product Quality: Our commitment to quality ensures reliable and reproducible results.

Competitive Pricing: We understand the importance of working within your research budget and offer competitive pricing.

Streamline your Supply Chain: Our products are produced in the U.S., are always in stock, and eliminate delays in your experiments.

Top-Tier Service & Support: Our teams of experts are dedicated to aid your project's success.

One-Year Guarantee: We stand behind the quality of our products, providing you with peace of mind as you focus on research.

Discover Faster: With Bio X Cell by your side, you can research with confidence and make immediate strides in your discoveries.

Product Categories



***InVivoMAB*[™] Antibodies**

High-quality, high-purity antibodies to drive success in your *in vivo* and *in vitro* applications.



***InVivoSIM*[™] Antibodies**

Study the biological effects of a drug with affordable research-grade biosimilar monoclonal antibodies.



***InVivoPlus*[™] Antibodies**

Exceed strict demands and rigorous standards required for *in vivo* work at any research organization.



RecombiMAB[™] Antibodies

Reduce immunogenicity in mouse models or humanized mouse models with our recombinant monoclonal antibodies.



Isotype Controls

Clearly identify antibody-specific signals with our isotype controls-delivered ready for *in vivo* administration.



***InVivoPure*[™] Dilution Buffer**

Preserve antibody stability with buffers specifically formulated and tested for *in vivo* applications.

in vivo Antibody Features

Exceptional Purity: Ultra-pure antibody solution without added proteins, preservatives, or chemicals. Each lot is QC tested for purity using SDS-PAGE.

Pathogen Free: All *InVivoPlus*[™] products are screened for an exhaustive panel of murine pathogens to help you adhere to IACUC and Animal Facility requirements.

Low Protein Aggregation: Antibody solution with low levels of protein aggregation. *InVivoPlus*[™] products are QC tested for aggregate level and guaranteed to be <5% of the total protein.

Ultra-low Endotoxin Levels: We QC test the endotoxin level for each lot. The endotoxin levels of *InVivoMAB*[™] products are < 2EU/mg and *InVivoPlus*[™] products are < 1EU/mg.

Advanced Binding Validation: We utilize a library of recombinant proteins to validate that each lot of applicable *InVivoPlus*[™] antibody binds strongly and specifically to its target antigen.

Matched Isotype Control Antibodies: We carry a wide selection of non-binding isotype control antibodies, taking the guesswork out of finding the correct control for your antibody.

Featured *InVivoMab*[™] and *InVivoPlus*[™] Antibodies

Our antibodies feature greater than 95% purity, ultra-low endotoxin levels, and are preservative, stabilizer, and carrier protein-free. They are used extensively in studying animal models of human diseases. With an ever-growing selection of functionally active, blocking, neutralizing, and agonistic antibodies, our products cover a wide range of research areas including **immuno-oncology**, **immunology**, **neuroscience** and more. Our products are cited in over 20,000 publications that detail a broad range of applications, including:

in vivo Applications

- Immune Checkpoint Blocking
- Mouse Cell Specific Depletion
- Cytokine Neutralization
- Agonistic / Cell Activation
- Receptor Blocking

in vitro Applications

- Functional Assays
- Neutralization Assays
- Cell-based Assays

Protein Detection

- WB, IP, IHC, IF, FC, ELISA

Antigen	Reactivity	Application	Clone	Catalog #	Isotype Control
CD3ε	Mouse	<i>in vivo</i> T cell depletion, <i>in vitro</i> T cell stimulation/activation, IF, FC, WB	145-2C11	BP0001	BP0091
CD4	Mouse	<i>in vivo</i> CD4+ T cell depletion	YTS 191	BE0119	BE0090
CD8α	Mouse	<i>in vivo</i> CD8+ T cell depletion, WB	2.43	BP0061	BP0090
CD8β	Mouse	<i>in vivo</i> CD8+ T cell depletion, <i>in vitro</i> CD8 blockade, IF	53-5.8	BE0223	BE0088
CD16/CD32	Mouse	<i>in vitro</i> Fc receptor blocking, <i>in vivo</i> Fc receptor blocking	2.4G2	BP0307	BP0090
CD19	Mouse	<i>in vivo</i> B cell depletion, <i>in vivo</i> CD19 neutralization, <i>in vitro</i> B cell negative selection, FC	1D3	BE0150	BE0089
CD20	Mouse	<i>in vivo</i> B cell depletion, WB	MB20-11	BP0356	BP0366
CD25	Mouse	<i>in vivo</i> regulatory T cell depletion, FC	PC-61.5.3	BP0012	BP0088
CD47	Human, Mouse, Rat	<i>in vivo</i> CD47 blockade, <i>in vitro</i> CD47 blocking, IF	MIAP410	BP0283	BP0083
CD71	Mouse	<i>in vivo</i> depletion of CD71+ cells	R17 217.1.3/TIB-219	BP0175	BP0089
CTLA-4	Mouse	<i>in vivo</i> CTLA-4 neutralization, WB	9D9	BP0164	BP0086
CTLA-4	Mouse	<i>in vivo</i> CTLA-4 neutralization, <i>in vitro</i> CTLA-4 neutralization, WB	9H10	BP0131	BP0087
IFNAR-1	Mouse	<i>in vivo</i> IFNAR-1 blockade, <i>in vitro</i> IFNAR-1 blockade, WB	MAR1-5A3	BP0241	BP0083
IFNγ	Mouse	<i>in vivo</i> IFNγ neutralization, <i>in vitro</i> IFNγ neutralization, ELISPOT, FC, WB	XMG1.2	BP0055	BP0088
IL-1β	Mouse, Rat	<i>in vivo</i> IL-1β neutralization, <i>in vitro</i> IL-1β neutralization, ELISA	B122	BE0246	BE0091
IL-4	Mouse	<i>in vivo</i> IL-4 neutralization, <i>in vitro</i> IL-4 neutralization, <i>in vivo</i> IL-4 receptor stimulation, FC, WB	11B11	BP0045	BP0088
IL-6	Mouse	<i>in vivo</i> IL-6 neutralization, <i>in vitro</i> IL-6 neutralization	MP5-20F3	BE0046	BE0088
IL-10R	Mouse	<i>in vivo</i> blocking of IL-10/IL-10R signaling, <i>in vitro</i> blocking of IL-10R signaling, FC, WB	1B1.3A	BP0050	BP0088
IL-12 p40	Mouse	<i>in vivo</i> IL-12p40 neutralization, p40 affinity chromatography, IP, ELISA, FC, WB	C17.8	BP0051	BP0089
IL-17A	Mouse	<i>in vivo</i> IL-17A neutralization	17F3	BP0173	BP0083
LAG-3	Mouse	<i>in vivo</i> LAG-3 neutralization, <i>in vitro</i> LAG-3 neutralization, FC, WB	C9B7W	BP0174	BP0088
LPAM-1	Mouse	<i>in vivo</i> Integrin α4β7 neutralization, FC	DATK32	BE0034	BE0089
Ly6G/Ly6C	Mouse	<i>in vivo</i> depletion of Gr-1+ myeloid cells, FC, IHC-P, IHC-F	RB6-8C5	BP0075	BP0090
MHC Class I	Human	Functional assays	W6/32	BP0079	BP0085
NK1.1	Mouse	<i>in vivo</i> NK cell depletion, FC	PK136	BP0036	BP0085
OX40	Mouse	<i>in vivo</i> OX40 activation, <i>in vitro</i> OX40 activation, WB	OX-86	BP0031	BP0088
PD-1	Mouse	<i>in vivo</i> blocking of PD-1/PD-L signaling	RMP1-14	BP0146	BP0089
PD-1	Mouse	<i>in vivo</i> blocking of PD-1/PD-L signaling, <i>in vitro</i> PD-1 neutralization, IHC-F, IF, WB, FC	29F.1A12 [™]	BP0273	BP0089
PD-L1	Mouse	<i>in vivo</i> PD-L1 blockade, IF, IHC-F, FC, WB	10F.9G2 [™]	BP0101	BP0090
TGF-β	Human, Mouse, Rat	<i>in vivo</i> TGFβ neutralization, <i>in vitro</i> TGFβ neutralization, WB	1D11.16.8	BP0057	BP0083
TNFα	Mouse	<i>in vivo</i> TNFα neutralization, <i>in vitro</i> TNFα neutralization, WB	XT3.11	BP0058	BP0088
VEGFR-2	Mouse	<i>in vivo</i> blocking of VEGF/VEGFR-2 signaling, <i>in vitro</i> blocking of VEGFR signaling, WB	DC101	BP0060	BP0088

For the full list of our products, visit us at bioxcell.com/in-vivo-antibodies.

Discover The *InVivoPlus*[™] Difference

The *InVivoPlus*[™] versions of our products are structurally and functionally identical to the *InVivoMab*[™] versions. The difference lies in a superior level of quality that will meet or exceed the strict demands and rigorous standards required for *in vivo* work at any research organization.

*Binding validation is performed for each lot for all applicable *InVivoPlus*[™] products. Non-applicable *InVivoPlus*[™] products include antibodies for which an antigen is not commercially available and non-reactive isotype control antibodies.

† Murine pathogen screening includes: Mycoplasma including *M. pulmonis*, *M. arginini*, *M. fermentans*, *M. hominis*, *M. hyorhinis*, *M. orale*, *M. pirum*, *M. salivarium*, *M. agassizii*, *M. cynos* and others. Murine norovirus (MNV), murine parvovirus (MPV) 1, 2, 3, 4, and 5, murine minute virus (MMV/MVM), murine hepatitis virus (MHV), murine reovirus (REO) type 1, 2, and 3, lymphocytic choriomeningitis virus (LCMV), lactate dehydrogenase elevating virus (LDV), murine rotavirus MRV/EDIM), Theiler's murine encephalomyelitis virus (TMEV), ectromelia virus (ECTRO), hantavirus (HANTA), polyoma virus (POLY), murine adenovirus (detects and differentiates MAD1 and MAD2), sendai virus (SEND), pneumonia virus of mice (PVM), murine cytomegalovirus (MCMV) 1 and 2 (MCMV), and K virus.

	<i>InVivoMab</i> [™]	<i>InVivoPlus</i> [™]
Binding Validation* Determined by immunoblot, flow cytometry, or ELISA		✓
Murine Pathogen Screening† Determined by ultrasensitive real-time PCR		✓
Antibody Aggregation Screening Determined by Size Exclusion Chromatography (SEC)		✓
Endotoxin Determined by LAL gel clotting assay	<2EU/mg	<1EU/mg
Purity Determined by SDS-PAGE	>95%	>95%
Preservative, stabilizer and carrier protein free	✓	✓
Formulated for <i>in vivo</i> use	✓	✓
Catalog number format	BE####	BP####

Featured *InVivoSIM*[™] Biosimilar Antibodies

The *InVivoSIM*[™] research-grade biosimilar antibodies have the same variable region sequences as the original therapeutic antibodies. Biosimilars make it possible to study the biological effects of a drug without the need to source an expensive pharmaceutical-grade therapeutic. They are an excellent choice for use as standard of care/therapeutic benchmarks in functional assays, pharmacokinetic assays, and *in vivo* studies in xenograft and humanized mouse models. [For the full list of our products, visit us at bioxcell.com/biosimilars.](https://bioxcell.com/biosimilars)

Product Name	Application	Catalog #	Isotype Control
<i>InVivoSIM</i> [™] anti-human C5 (Eculizumab Biosimilar)	Inhibition of the activation of C5, Functional assays, ELISA, IP	SIM0011	CP147
<i>InVivoSIM</i> [™] anti-human CD20 (Rituximab Biosimilar)	FC, ELISA, WB	SIM0008	BP0297
<i>InVivoSIM</i> [™] anti-human CTLA-4 (Ipilimumab Biosimilar)	<i>in vivo</i> CTLA-4 neutralization, FC, ELISA, WB	SIM0004	BP0297
<i>InVivoSIM</i> [™] anti-human EGFR (Cetuximab Biosimilar)	<i>in vitro</i> EGFR blockade, EGFR blockade, ELISA, FC	SIM0002	BP0297
<i>InVivoSIM</i> [™] anti-human HER2 (Trastuzumab Biosimilar)	FC, ELISA, IHC, WB	SIM0005	BP0297
<i>InVivoSIM</i> [™] anti-human HER2 (Pertuzumab Biosimilar)	Reported Applications, Functional assays, ELISA, Neutralization	SIM0019	BE0297
<i>InVivoSIM</i> [™] anti-human IgE (Omalizumab Biosimilar)	IgE neutralization, FC, ELISA	SIM0016	BP0297
<i>InVivoSIM</i> [™] anti-human IL-12 p40 (Ustekinumab Biosimilar)	Functional assays	SIM0020	BP0297
<i>InVivoSIM</i> [™] anti-human IL-17A (Secukinumab Biosimilar)	Functional assays, ELISA, IHC, IHC, FC	SIM0013	BP0297
<i>InVivoSIM</i> [™] anti-human IL-6R (Tocilizumab Biosimilar)	Functional assays, ELISA	SIM0014	BP0297
<i>InVivoSIM</i> [™] anti-human LAG-3 (Relatlimab Biosimilar)	ELISA	SIM0015	CP147
<i>InVivoSIM</i> [™] anti-human NGF (Tanezumab Biosimilar)	NGF neutralization, ELISA	SIM0017	BP0301
<i>InVivoSIM</i> [™] anti-human PD-1 (Nivolumab Biosimilar)	Blocking of PD-1/PD-L signaling, FC, IHC, WB	SIM0003	CP147
<i>InVivoSIM</i> [™] anti-human PD-1 (Pembrolizumab Biosimilar)	<i>in vivo</i> Blocking of PD-1/PD-L signaling, Functional assays	SIM0010	CP147
<i>InVivoSIM</i> [™] anti-human PD-L1 (Atezolizumab Biosimilar)	<i>in vitro</i> PD-L1 blockade, FC, WB	SIM0009	BP0297
<i>InVivoSIM</i> [™] anti-human PD-L1 (Avelumab Biosimilar)	<i>in vitro</i> PD-L1 blockade, Functional assays, ELISA	SIM0021	CP169
<i>InVivoSIM</i> [™] anti-human TNF α (Adalimumab Biosimilar)	TNF α neutralization, FC, ELISA, IHC, IP, IHC, WB	SIM0001	BP0297
<i>InVivoSIM</i> [™] anti-human TNF α (Infliximab Biosimilar)	<i>in vitro</i> TNF α neutralization, FC, ELISA, WB	SIM0006	BP0297
<i>InVivoSIM</i> [™] anti-human VEGF (Bevacizumab Biosimilar)	VEGF neutralization, FC, ELISA, IP, WB	SIM0007	BP0297
<i>InVivoSIM</i> [™] anti-human VEGFR-2 (Ramucirumab Biosimilar)	Functional assays, IHC, FC	SIM0012	BP0297

Featured RecombiMab[™] Recombinant Antibodies

RecombiMab[™] antibodies use mouse or human IgG constant regions instead of the typical rat or hamster IgG constant regions, resulting in reduced immunogenicity in mouse and humanized mouse models. In some syngeneic mouse tumor models, particularly in BALB/c mice, repeated rat or hamster IgG administration can result in complications from hypersensitivity reactions. Using the RecombiMab[™] antibody instead of wild-type clones may reduce antibody immunogenicity and alleviate hypersensitivity reactions. RecombiMab[™] antibodies also overcome the limitations of traditionally manufactured antibodies by ensuring high lot-to-lot consistency and data reproducibility across experiments. [For the full list of our products, visit us at bioxcell.com/recombimab.](https://bioxcell.com/recombimab)

Antigen	Reactivity	Host / Isotype	Mutation(s)	Application(s)	Clone	Catalog #	Isotype Control
CD4	Mouse	Mouse, IgG2b, κ		<i>in vivo</i> CD4+ T cell depletion, FC, WB	GK1.5-CP127	CP127	BP0086
CD8 α	Mouse	Mouse, IgG2b, κ		<i>in vivo</i> CD8+ T cell depletion, WB	2.43-CP128	CP128	BP0086
CD8 α	Mouse	Mouse, IgG2a, κ		<i>in vivo</i> CD8+ T cell depletion, WB	YTS 169.4-CP134	CP134	BP0085
CD16/CD32	Mouse	Mouse, IgG2a		<i>in vitro</i> Fc receptor blocking, <i>in vivo</i> Fc receptor blocking	2.4G2-CP025	CP025	BP0085
CD40	Mouse	Mouse, IgG2a, κ		<i>in vivo</i> CD40 activation, <i>in vitro</i> B cell stimulation/activation	FGK4.5-CP133	CP133	BP0085
CD71 (TfR1)	Mouse	Mouse, IgG2a, κ		<i>in vivo</i> depletion of CD71+ cells	R17 2171.3-CP130	CP130	BP0085
CSF1R (CD115)	Mouse	Mouse, IgG2a, κ		<i>in vivo</i> macrophage & monocyte depletion, <i>in vitro</i> CSF1R neutralization, FC, WB	AFS98-CP131	CP131	BP0085
CTLA-4 (CD152)	Mouse	Mouse, IgG1, κ		<i>in vivo</i> CTLA-4 neutralization, <i>in vitro</i> CTLA-4 neutralization, WB	9H10-CP146	CP146	BP0083
CTLA-4 (CD152)	Mouse	Mouse, IgG1		<i>in vivo</i> CTLA-4 neutralization, WB	9D9-CP006	CP006	BP0083
LAG-3	Mouse	Mouse, IgG2a		<i>in vivo</i> LAG-3 neutralization, <i>in vitro</i> LAG-3 neutralization, FC, WB	C9B7W-CP014	CP014	BP0085
OX40 (CD134)	Mouse	Mouse, IgG2a		<i>in vivo</i> OX40 activation*, <i>in vitro</i> OX40 activation*, Western blot*	OX86-CP017	CP017	BP0085
PD-1 (CD279)	Mouse	Mouse, IgG2a, κ	D265A	<i>in vivo</i> blocking of PD-1/PD-L signaling	RMP1-14-CP151	CP151	CP150
PD-1 (CD279)	Mouse	Mouse, IgG2a, κ	L234A, L235A, P329G (LALA-PG)	<i>in vivo</i> blocking of PD-1/PD-L signaling	RMP1-14-CP153	CP153	CP150
PD-1 (CD279)	Mouse	Mouse, IgG2a, κ		<i>in vivo</i> blocking of PD-1/PD-L signaling	RMP1-14-CP157	CP157	BP0085
PD-1 (CD279)	Mouse	Mouse, IgG1, κ		<i>in vivo</i> blocking of PD-1/PD-L signaling, <i>in vitro</i> PD-1 neutralization, IHC-F, FC, WB	29F.1A12-CP159	CP159	BP0083
PD-1 (CD279)	Mouse	Mouse, IgG1, κ		<i>in vivo</i> blocking of PD-1/PD-L signaling	RMP1-14-CP162	CP162	BP0083
PD-1 (CD279)	Mouse	Mouse, IgG1		<i>in vivo</i> blocking of PD-1/PD-L signaling, <i>in vitro</i> PD-1 neutralization, IHC-F, IF, WB, FC	29F.1A12-CP004	CP004	BP0083
PD-L1 (B7-H1)	Mouse	Mouse, IgG1, κ		<i>in vivo</i> PD-L1 blockade, IF, IHC-F, FC	10F.9G2-CP168	CP168	BP0083
VEGFR-2	Mouse	Mouse, IgG2a, κ		<i>in vivo</i> blocking of VEGF/VEGFR-2 signaling, <i>in vitro</i> blocking of VEGFR signaling, WB	DC101-CP132	CP132	BP0085

Contract Services

Enlist the Experts at Bio X Cell for Antibody Production Services

At Bio X Cell, we know that growing hybridoma cells and purifying antibodies in large quantities with low endotoxin levels can be technically challenging and extremely time-consuming. Our optimized fermentation process allows us to scale up antibody production from milligrams to grams much faster than traditional methods. A technical project manager will guide you through the process, providing recommendations and updates along the way. We know that research moves at a rapid and precise pace, and we aim to partner with you to achieve your goals faster.

To discuss your contract services projects, please reach us at contractservices@bioxcell.com

- 25+ years of expertise
- Customizable services
- Low endotoxin levels
- Fast turnaround time
- Ultra-high purity
- Large scale production

Antibody Production Services: We know that growing hybridoma cells and purifying antibodies can be technically challenging and time-consuming. Enlist our experts to do this for you so you can focus on your research.

Antibody Conjugation Services: Gain access to a wide selection of antibodies and labels. Quickly get high-quality conjugated antibodies ready and specific for your application - saving you time and cost.

Transient Expression Services: Large scale transient recombinant antibody expression services with fast turnaround time to support your *in vivo* and *in vitro* studies.

Cell Freezing and Storage: Optimal storage that you can trust. We have optimized processes to create master and working cell banks to ensure your cells are safe and readily available for scale-up.

Murine Pathogen Testing: We use ultra-sensitive qPCR to screen for a multitude of murine pathogens, ensuring that the antibodies you are using for *in vivo* research will meet all requirements of your research facility.

Mycoplasma Testing: Don't let mycoplasma jeopardize your valuable cell line. We offer services to reliably detect mycoplasma and restore your cells to a usable state to keep your projects on track.



Impacting the World Beyond Antibodies

Building a Green and Sustainable Business

At Bio X Cell we are committed to advancing our mission with minimal environmental impact. Our employee advisory committee convenes regularly to explore the ways in which we can go green as a company. In 2024, we transitioned to recyclable and biodegradable plastic packaging for domestic outbound shipments. This eco-friendly alternative retains optimal thermal properties during shipping, saves thousands of pounds of waste, and cuts down on our carbon emissions annually. In 2023, we installed both ground and rooftop solar arrays at our headquarters. By harnessing sunlight and transforming it into a sustainable energy source, we are significantly reducing our carbon footprint and contributing to a cleaner environment.

Our Philanthropy

The Bio X Cell team is united by a shared passion for giving back to our community and enhancing the lives of those around us. In 2019, we launched the Bio X Cell Fund, a philanthropic initiative driven by our employees and managed in partnership with the New Hampshire Charitable Foundation. The Bio X Cell Fund supports organizations that address mental and physical health issues, environmental issues, education, arts and culture, and humanitarian aid, including disaster relief. Over the past five years, we have raised over \$1.24 million for various foundations, creating a positive community impact that enriches both our professional and personal lives.





For the full list of our products, scan the QR code or visit us at bioxcell.com

UNITED STATES
BIO X CELL HEADQUARTERS
39 Labombard Rd, Lebanon, NH 03766
Tel: 1-866-787-3444
customerservice@bioxcell.com
sales@bioxcell.com
www.bioxcell.com

EUROPE
BIO X CELL EUROPEAN SUPPLY CENTER
Tel: +49 6171 8669416
Tel (mobile): +49 17681170500
euro-orders@bioxcell.com
www.bioxcell.com

CHINA
BIO X CELL CHINA
sales@bioxcell.com.cn
www.bioxcell.com.cn

KOREA
BIO X CELL KOREA
Tel: 02-2013-8240
bioxcell@bioxcell.co.kr
www.bioxcell.co.kr

INTERNATIONAL DISTRIBUTORS

For a full list of international distributors, visit: bioxcell.com/our-distributors/

For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Bio X Cell, Bio X Cell Logo and all other trademarks are property of Bio X Cell.
© 2024 Bio X Cell

